

# See the difference **reditest**® can make.

Equip your drug test program with reliable and rapid on-site screening devices.

To place an order, call: **1-877-444-0049**. *Can't find what your looking for? Just ask us!*

With Reditest, we offer you more than the highest quality, most cost-effective devices in the industry.

We offer complete drug testing confidence.

## Reditest® Panel-Dip URINE SUBSTANCE ABUSE SCREENING DEVICES

PANEL(S)	PANEL-DIP CONFIGURATIONS — 25 per box	PART NUMBER
1	AMPHETAMINES (AMP)	01 102 0018
1	METHAMPHETAMINES (M-AMP)	01 102 0002
1	BARBITURATES (BAR)	01 102 0019
1	BUPRENORPHINE (BUP)	01 102 1955
1	BENZODIAZEPINES (BZO)	01 102 0022
1	COCAINE (COC)	01 102 0001
1	METHADONE (MTD)	01 102 0020
1	MDMA (Ecstasy)	01 102 0036
1	OPIATES (300) (OPI)	01 102 0003
1	OXYCODONE (OXY)	01 102 0037
1	TRICYCLIC ANTIDEPRESSANTS (TCA)	01 102 0023
1	MARIJUANA (THC)	01 102 0004
1	PHENCYCLIDINE (PCP)	01 102 0021
2	COC/M-AMP	01 102 0007
2	COC/OPI(300)	01 102 0005
2	COC/THC	01 102 0006
2	M-AMP/THC	01 102 0008
2	M-AMP/OPI(300)	01 102 0030
3	COC/M-AMP/THC	01 102 0009
3	COC/THC/OPI(300)	01 102 0010
3	M-AMP/THC/OPI(300)	01 102 0011
3	M-AMP/COC/OPI(300)	01 102 0014
4	COC/M-AMP/THC/OPI(300)	01 102 0012
4	COC/AMP/THC/OPI(300)	01 102 0032
5	COC/M-AMP/THC/OPI(300)/BZO	01 102 0015
5	COC/M-AMP/THC/OPI(300)/PCP	01 102 0013
5	COC/AMP/THC/OPI(300)/PCP	01 102 0033
5	COC/AMP/M-AMP/THC/OPI(300)	01 102 0034
5	COC/AMP/THC/OPI2(2000)/PCP	01 102 0047
6	COC/M-AMP/THC/OPI(300)/PCP/BZO	01 102 0016
6	COC/M-AMP/THC/MTD/OPI(300)/BZO	01 102 0017
6	COC/M-AMP/THC/OPI(300)/BAR/BZO	01 102 0024
6	COC/M-AMP/THC/OPI(300)/OXY/BZO	01 102 0119
7	COC/AMP/THC/OPI(300)/PCP/BZO/TCA	01 102 0035
8	COC/AMP/M-AMP/THC/OPI(300)/OXY/BZO/MDMA	01 102 0169
10	COC/AMP/M-AMP/THC/MTD/OPI(300)/PCP/BAR/BZO/TCA	01 102 0025
10	COC/M-AMP/THC/MTD/MDMA/OPI(300)/OXY/PCP/BAR/BZO	01 102 0138
12	COC/AMP/M-AMP/THC/MTD/MDMA/OPI(300)/OXY/PPX/PCP/BAR/BZO	01 102 0141

## Reditest® RediCup® URINE SUBSTANCE ABUSE SCREENING DEVICES

PANEL(S)	REDICUP® CONFIGURATIONS — 25 per box	PART NUMBER
4	M-AMP/COC/THC/OPI(300)	01 102 0026
5	M-AMP/COC/THC/OPI(300)/BZO	01 102 0027
5	M-AMP/COC/THC/OPI(300)/PCP	01 102 0028
5	M-AMP/COC/THC/OPI(300)/AMP	01 102 0121
6	M-AMP/COC/THC/AMP/OPI2(2000)/BZO	01 102 0135
6	M-AMP/COC/THC/OPI(300)/PCP/BZO	01 102 0029
10	AMP/COC/THC/M-AMP/OPI2(2000)/PCP/BZO/TCA/BAR/MTD	01 102 0058
10	M-AMP/COC/THC/BZO/TCA/BAR/PCP/AMP/OPI(300)/MTD	01 102 0059
10	M-AMP/COC/THC/BZO/MTD/BAR/MDMA/OPI(300)/PCP/OXY	01 102 0137

## iScreen® iCup® URINE SUBSTANCE ABUSE SCREENING DEVICES

PANEL(S)	ICUP A.D CONFIGURATIONS (with Specimen Validity Test) — 25 per box	PART NUMBER
4	COC/THC/OPI(2000)/M-AMP—CLIA Waived	01 102 2032
4	COC(150)/THC/AMP/M-AMP(500)	01 102 2033
5	COC/THC/OPI(2000)/AMP/M-AMP—CLIA Waived	01 102 2034
5	COC/THC/MOP(300)/AMP/M-AMP	01 102 2021
6	COC/THC/OPI(2000)/AMP/M-AMP/BZO	01 102 2022
6	COC/THC/OPI(2000)/AMP(300)/OXY/MDMA	01 102 2037
8	COC/THC/OPI(2000)/AMP/M-AMP/PCP/BZO/BAR	01 102 2038
8	COC/THC/MOP(300)/AMP/M-AMP/BZO/OXY/PCP	01 102 2069
9	COC/THC/OPI(2000)/AMP/M-AMP/PCP/BZO/BAR/MTD	01 102 2039
10	COC/THC/OPI(2000)/AMP/M-AMP/PCP/BZO/BAR/MTD/TCA	01 102 2129
10	COC/THC/OPI(2000)/AMP/M-AMP/BZO/BAR/OXY/MDMA/PPX	01 102 2074
12	COC/THC/OPI(2000)/AMP/M-AMP/PCP/BZO/BAR/MTD/TCA/OXY/PPX	01 102 2027
PANEL	ICUP ALL INCLUSIVE CUP CONFIGURATION	PART NUMBER
10	COC/THC/OPI(2000)/AMP/M-AMP/BZO/BAR/OXY/MDMA/PPX	01 102 2020

## E-Z Split Key Cup II URINE SUBSTANCE ABUSE SCREENING DEVICES

PANEL	EZ SPLIT KEY CUP II CONFIGURATION (with Specimen Validity Test) — 25 per box	PART NUMBER
5	COC/THC/OPI(2000)/AMP/PCP—CLIA Waived	01 102 1974
PANEL(S)	EZ SPLIT KEY CUP II CONFIGURATIONS — 25 per box	PART NUMBER
6	COC/THC/OPI(2000)/AMP/M-AMP/BZO	01 102 1984
10	COC/THC/OPI(2000)/AMP/M-AMP/PCP/BZO/BAR/MTD/MDMA	01 102 1985
12	COC/THC/MOP(300)/AMP/M-AMP/BZO/BAR/MTD/MDMA/OXY/BUP/PPX	01 102 2096

## Reditest® 6 Cassette CLIA WAIVED URINE SUBSTANCE ABUSE SCREENING DEVICE

PANEL	REDITEST® 6 CASSETTE CONFIGURATION — 25 per box	PART NUMBER
6	AMP/COC/M-AMP/OPI2(2000)PCP/THC—CLIA Waived	01 102 0166

## Reditest® On-Site Oral ORAL FLUID SUBSTANCE ABUSE SCREENING DEVICE

PANEL	REDITEST® ON-SITE ORAL CONFIGURATION — 25 per box	PART NUMBER
6	AMP/COC/M-AMP/OPI2(2000)/PCP/THC	01 102 0127

## iScreen® OFD ORAL FLUID SUBSTANCE ABUSE SCREENING DEVICES

PANEL(S)	ISCREEN® OFD CONFIGURATIONS — 25 per box	PART NUMBER
5	COC/THC/OPI(2000)/M-AMP/AMP	01 102 2024
6	COC/THC/OPI(2000)/M-AMP/AMP/PCP	01 102 2025

## OrAlert™ ORAL FLUID SUBSTANCE ABUSE SCREENING DEVICES

PANEL(S)	ORALERT CONFIGURATIONS — 25 per box	PART NUMBER
6	AMP/COC/M-AMP/OPI(2000)/PCP/THC	01 102 1960
6	AMP/COC/M-AMP/OPI(2000)/BZO/THC	01 102 2083

## Reditest® Alcohol-Saliva Test ORAL FLUID ALCOHOL SCREENING DEVICE

DEVICE	REDITEST® ALCOHOL-SALIVA TEST STRIP — 25 per box	PART NUMBER
TEST STRIP	DETECTS 0.02%/ 0.08%/ AND 0.30% SAC	01 362 0001

## Reditest® Alcohol-Breath Test BREATH ALCOHOL SCREENING DEVICE

DEVICE	REDITEST® ALCOHOL-BREATH TEST DEVICE — 25 per box	PART NUMBER
TEST STRIP	DETECTS 0.02%/ 0.08%/ AND 0.30% SAC	01 215 0004

## ALCO-SCREEN® ORAL FLUID ALCOHOL SCREENING DEVICES

DEVICE	ALCO-SCREEN® // ALCO-SCREEN® 02 — 25 per box	PART NUMBER
TEST STRIP	DETECTS 0.02%/ 0.08%/ AND 0.30% SAC	01 094 0055
TEST STRIP	DETECTS 0.02%—DOT APPROVED	01 094 0056

## Reditest® Smoke Cassette URINE COTININE SCREENING DEVICES

DEVICE	REDITEST® SMOKE CASSETTE — 25 per box	PART NUMBER
CASSETTE	COTININE (200 NG/ML)—DETECTS NICOTINE METABOLITE	01 102 0140

## Adulteration Strips SPECIMEN VALIDITY TEST (SVT)

DEVICE	ADULTERATION STRIPS — 25 strips per canister	PART NUMBER
TEST STRIP	1-STEP VALIDITY TEST. SCREENS 6 PARAMETERS FOR SPECIMEN VALIDITY, DILUTION & TAMPERING. <i>Creatinine, Nitrite, Glutaraldehyde, pH, Specific Gravity, and Oxidants/Pyridinium Chlorochromate</i>	01 102 1910

# screening device cut-off levels

The Reditest devices detect a wide range of drugs and drug metabolites in urine or oral fluids at the following cut-off concentrations:

## URINE SCREENING DEVICE CUT-OFF LEVELS

DRUG	CALIBRATOR	CUT-OFF LEVEL
Amphetamine (AMP 1000)	d-Amphetamine	*1,000 ng/mL
Amphetamine (AMP 300)	d-Amphetamine	300 ng/mL <sup>1</sup>
Barbiturates (BAR)	Secobarbital	300 ng/mL
Benzodiazepines (BZO)	Oxazepam	300 ng/mL
Buprenorphine (BUP)	Buprenorphine	10 ng/mL
Cocaine (COC 300)	Benzoyllecgonine	*300 ng/mL
Cocaine (COC 150)	Benzoyllecgonine	150 ng/mL <sup>2</sup>
Methylenedioxymethamphetamine (MDMA)	D,L Methylenedioxymethamphetamine	500 ng/mL
Marijuana (THC)	11-nor- $\Delta^9$ -THC-9 COOH	*50 ng/mL
Methadone (MTD)	Methadone	300 ng/mL
Methamphetamine (M-AMP)	d-Methamphetamine	1000 ng/mL
Methamphetamine (M-AMP 500)	d-Methamphetamine	500 ng/mL <sup>2</sup>
Opiates (OPI 300 / MOP 300)	Morphine	300 ng/mL
Opiates (OPI2 2000 / OPI 2000)	Morphine	2000 ng/mL
Oxycodone (OXY)	Oxycodone	100 ng/mL
Phencyclidine (PCP)	Phencyclidine	*25 ng/mL
Propoxyphene (PPX)	Propoxyphene	300 ng/mL
Tricyclic Antidepressants (TCA)	Nortriptyline	1000 ng/mL

\*Meets SAMHSA cut-off levels.

## SALIVA SCREENING DEVICE CUT-OFF LEVELS

DRUG	CALIBRATOR	CUT-OFF LEVEL
Amphetamine (AMP)	d-Amphetamine	50 ng/mL
Methamphetamine (M-AMP)	d-Methamphetamine	50 ng/mL
Benzodiazepines (BZO)	Oxazepam	10 ng/mL
Cocaine (COC)	Benzoyllecgonine	20 ng/mL
Opiates (OPI)	Morphine	40 ng/mL
Marijuana (THC)	11-nor- $\Delta^9$ -THC-9 COOH	12 ng/mL
Marijuana (THC)	$\Delta^9$ -THC	100 ng/mL <sup>3</sup>
Phencyclidine (PCP)	Phencyclidine	10 ng/mL

1. iScreen® iCup® P/N: 01 102 2037, AMP detection at 300 ng/mL.

2. iScreen® iCup® P/N: 01 102 2033, COC detection at 150 ng/mL, AMP detection at 500 ng/mL.

3. OrAlert THC detection. Includes improved sensitivity of parent THC.

RediCup®, Panel-Dip, 6 Cassette devices, EZ-Split Key Cup, iCup® and Alcohol-Breath Test are for professional in vitro diagnostic use only. On-Site Oral, Alcohol-Saliva Strips, iScreen™ OFD, OrAlert™ and ALCO-SCREEN are for forensic use only. Reditest® Smoke Cassette is for determination of smoking status only. ALCO-SCREEN O2 is for DOT use. The assays provide only a preliminary test result. Any positive result obtained with these drug screens are presumptive and should be confirmed by an alternate method such as GC/MS.

11 000 0135 REV15